

Hypersil GOLD HPLC and UHPLC Columns

Outstanding peak shape for your separations

Thermo Scientific™ Hypersil GOLD™ columns are an ideal choice for new method development, delivering outstanding peak shape and increased sensitivity for your separations.

- **Solutions for all Separation Needs**
Particle sizes and column designs optimized for all analyses
- **Improved Sensitivity**
Excellent peak symmetry
- **Enhanced Resolution**
Separation of complex samples
- **pH Stability**
Robust assays at extended pH




The Hypersil GOLD web page contains the latest news, applications and downloads for the Hypersil GOLD HPLC column range. Visit www.thermoscientific.com/hypersilgold


Hypersil GOLD

Endcapped, ultra-pure, silica-based columns with exceptional peak shape and resolution for HPLC and LC-MS

- Significant reduction in peak tailing while retaining C18 selectivity
- Excellent resolution, efficiency and sensitivity
- Confidence in the accuracy and quality of analytical data



Seven β-blockers in 1 minute



Hypersil GOLD, 1.9µm, 20 x 2.1mm

Mobile Phase A:	H ₂ O+0.1%formic acid
Mobile Phase B:	MeCN+0.1%formic acid
Gradient:	15 to 100% B in 1min
Temperature:	30°C
Flow Rate:	0.5mL/min
Detection:	+ESI
Analytes:	1. Atenolol 2. Nadolol 3. Pindolol 4. Timolol 5. Metoprolol 6. Oxprenolol 7. Propranolol

Hypersil GOLD

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID	
1.9	UHPLC Column	20	-	25002-022130	-	-	-	
		30	25002-031030	25002-032130	-	-	-	
		50	25002-051030	25002-052130	25002-053030	-	25002-054630	
		100	25002-101030	25002-102130	25002-103030	-	-	
		150	-	25002-152130	-	-	-	
		200	-	25002-202130	-	-	-	
3	Drop-in Guard (4/pk)	10	25003-011001	25003-012101	25003-013001	25003-014001	25003-014001	
		30	-	25003-032130	-	-	25003-034630	
	HPLC Column	50	-	25003-052130	25003-053030	25003-054030	25003-054630	
		100	25003-101030	25003-102130	25003-103030	25003-104030	25003-104630	
		150	25003-151030	25003-152130	25003-153030	25003-154030	25003-154630	
		200	-	25003-202130	-	-	-	
5	Drop-in Guard (4/pk)	10	-	25005-012101	25005-013001	25005-014001	25005-014001	
		30	-	25005-032130	25005-033030	-	25005-034630	
	HPLC Column	50	-	25005-052130	25005-053030	-	25005-054630	
		100	-	25005-102130	25005-103030	-	25005-104630	
		150	-	25005-152130	25005-153030	25005-154030	25005-154630	
		200	-	25005-202130	25005-203030	25005-204030	25005-204630	
		UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00

Hypersil GOLD Preparative



Particle Size (µm)	Format	Length (mm)	10mm ID	20mm ID	21mm ID	30mm ID	50mm ID
5	Preparative Guard Cartridge (3/pk)	10	25005-019023	25005-019223	-	-	-
	Preparative HPLC Column	50	25005-059070	-	-	25005-059370	25005-059570
		100	25005-109070	-	25005-109270	25005-109370	25005-109570
		150	25005-159070	-	25005-159270	25005-159370	25005-159570
		250	25005-259070	-	25005-259270	25005-259370	25005-259570
12	Preparative Guard Cartridge (3/pk)	10	25012-019023	25012-019223	-	-	-
	Preparative HPLC Column	50	25012-059070	-	25012-059270	25012-059370	25012-059570
		100	25012-109070	-	25012-109270	25012-109370	25012-109570
		150	-	-	25012-159270	25012-159370	25012-159570
		250	25012-259070	-	25012-259270	25012-259370	25012-259570
	Preparative Guard Cartridge Holder	10	C-1000	854-00	-	-	-



Hypersil GOLD C8

Recommended for analytes with medium hydrophobicity or when a less hydrophobic phase is required to obtain optimum retention

- Similar selectivity to C18 columns but with reduced retention
- Lower hydrophobicity, allowing compounds to elute quicker
- Faster separations
- Excellent peak shape
- High efficiency
- Outstanding sensitivity



Hypersil GOLD C8 5μm, 150 x 4.6mm

Mobile Phase:	MeOH
Temperature:	25°C
Flow Rate:	1.5mL/min
Detection:	UV, 450nm
Analytes:	1. Lutein 2. Lycopene 3. β-Carotene


Hypersil GOLD C8

Particle Size (μm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	20	-	25202-022130	-	-	-
		30	-	25202-032130	-	-	-
		50	-	25202-052130	25202-053030	-	25202-054630
		100	25202-101030	25202-102130	25202-103030	-	-
		150	-	25202-152130	-	-	-
3	Drop-in Guard (4/pk)	10	25203-011001	25203-012101	25203-013001	25203-014001	25203-014001
	HPLC Column	50	-	25203-052130	25203-053030	25203-054030	25203-054630
		100	25203-101030	25203-102130	25203-103030	25203-104030	25203-104630
		150	-	25203-152130	25203-153030	25203-154030	25203-154630
5	Drop-in Guard (4/pk)	10	-	25205-012101	25205-013001	25205-014001	25205-014001
	HPLC Column	50	-	25205-052130	25205-053030	-	25205-054630
		100	-	25205-102130	25205-103030	-	25205-104630
		150	-	25205-152130	25205-153030	25205-154030	25205-154630
		250	-	25205-252130	25205-253030	25205-254030	25205-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00


Hypersil GOLD C4

Lower hydrophobicity than C18 or C8 recommended for very hydrophobic analytes

- Lower hydrophobicity
- Faster separations
- Excellent peak shape
- High efficiency
- Outstanding sensitivity



Fatty acids



Hypersil GOLD C4 1.9µm, 100 x 2.1mm

Mobile Phase:	H ₂ O / MeCN (20:80)
Temperature:	30°C
Flow Rate:	0.55mL/min
Injection Volume:	1µL
Detection:	200 nm
Analytes:	1. Linolenic acid 2. Linoleic acid


Hypersil GOLD C4

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	25502-052130	-	-
		100	-	25502-102130	-	-
		150	-	25502-152130	-	-
3	Drop-in Guard (4/pk)	10	25503-011001	25503-012101	25503-013001	25503-014001
	HPLC Column	50	-	25503-052130	-	-
		100	-	25503-102130	25503-103030	25503-104630
150		25503-151030	25503-152130	25503-153030	25503-154630	
5	Drop-in Guard (4/pk)	10	-	25505-012101	25505-013001	25505-014001
		HPLC Column	50	-	25505-052130	-
	HPLC Column	100	-	25505-102130	25505-103030	25505-104630
		150	-	25505-152130	-	25505-154630
		250	-	25505-252130	25505-253030	25505-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00


Hypersil GOLD aQ

Controlled interaction mechanism by which polar analytes can be retained and resolved

- Polar endcapped C18 phase for alternative selectivity
- Retention and resolution of polar analytes
- Stable in 100% aqueous mobile phases



Water soluble vitamins



Hypersil GOLD aQ 5µm, 150 x 4.6mm

Mobile Phase A:	50 mM KH ₂ PO ₄ , pH 3.5
Mobile Phase B:	MeOH
Gradient:	0 – 100% B in 15 min
Flow Rate:	1mL/min
Detection:	UV, 205nm
Analytes:	<ol style="list-style-type: none"> 1. Vitamin B1 (thiamine) 2. Vitamin B6 (pyridoxine) 3. Vitamin B3 (nicotinamide) 4. Vitamin B5 (pantothenic acid) 5. Folic Acid 6. Vitamin B12 (cyanocobalamin) 7. Vitamin H (biotin) 8. Vitamin B2 (riboflavin)

Hypersil GOLD aQ

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	20	-	25302-022130	-	-	-
		30	-	25302-032130	-	-	-
		50	25302-051030	25302-052130	25302-053030	-	25302-054630
		100	25302-101030	25302-102130	25302-103030	-	-
		150	-	25302-152130	-	-	-
		200	-	25302-202130	-	-	-
3	Drop-in Guard (4/pk)	10	25303-011001	25303-012101	25303-013001	25303-014001	25303-014001
		30	-	25303-032130	25303-033030	-	-
	HPLC Column	50	-	25303-052130	25303-053030	25303-054030	25303-054630
		100	25303-101030	25303-102130	25303-103030	25303-104030	25303-104630
		150	25303-151030	25303-152130	25303-153030	25303-154030	25303-154630
		200	-	25303-202130	-	-	-
5	Drop-in Guard (4/pk)	10	-	25305-012101	25305-013001	25305-014001	25305-014001
		30	-	25305-032130	25305-033030	-	-
	HPLC Column	50	-	25305-052130	25305-053030	-	25305-054630
		100	-	25305-102130	25305-103030	-	25305-104630
		150	-	25305-152130	25305-153030	-	25305-154630
		250	-	25305-252130	-	25305-254030	25305-254630
UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00	


Hypersil GOLD aQ Preparative

Particle Size (µm)	Format	Length (mm)	10mm ID	20mm ID	21mm ID	30mm ID	50mm ID
5	Preparative Guard Cartridge (3/pk)	10	25305-019023	25305-019223	-	-	-
		50	25305-059070	-	25305-059270	25305-059370	25305-059570
	Preparative HPLC Column	100	25305-109070	-	25305-109270	25305-109370	25305-109570
		150	25305-159070	-	25305-159270	25305-159370	25305-159570
		250	25305-259070	-	25305-259270	25305-259370	25305-259570


Hypersil GOLD PFP

Introduction of a fluorine group into the stationary phase causes significant changes in solute-stationary phase interaction

- Fluorine atoms around the phenyl ring enhance pi-pi interactions with aromatic molecules
- Extra retention for halogenated species
- Selectivity for non-halogenated polar compounds



Polyphenols



Hypersil GOLD PFP 1.9μm, 50 x 2.1mm

Mobile Phase:	0.1% Acetic Acid
Temperature:	25°C
Flow Rate:	0.5mL/min
Injection Volume:	0.5μL
Detection:	UV, 280nm
Analytes:	1. Pyrogallol 2. Hydroquinone 3. Resorcinol 4. Phenol

Hypersil GOLD PFP

Particle Size (μm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	20	-	25402-022130	-	-	-
		30	-	25402-032130	-	-	-
		50	25402-051030	25402-052130	25402-053030	-	25402-054630
		100	25402-101030	25402-102130	25402-103030	-	-
		150	-	25402-152130	-	-	-
		200	-	25402-202130	-	-	-
3	Drop-in Guard (4/pk)	10	25403-011001	25403-012101	25403-013001	25403-014001	25403-014001
	HPLC Column	30	-	25403-032130	25403-033030	-	25403-034630
		50	-	25403-052130	25403-053030	25403-054030	25403-054630
		100	25403-101030	25403-102130	25403-103030	25403-104030	25403-104630
		150	25403-151030	25403-152130	25403-153030	-	25403-154630
5	Drop-in Guard (4/pk)	10	-	25405-012101	25405-013001	25405-014001	25405-014001
	HPLC Column	50	-	25405-052130	25405-053030	-	-
		100	-	25405-102130	25405-103030	-	25405-104630
		150	-	25405-152130	25405-153030	25405-154030	25405-154630
		250	-	25405-252130	-	25405-254030	25405-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00


Hypersil GOLD PFP Preparative

Particle Size (μm)	Format	Length (mm)	10mm ID	20mm ID	21mm ID	30mm ID	50mm ID
5	Preparative Guard Cartridge (3/pk)	10	25405-019023	25405-019223	-	-	-
	Preparative HPLC Column	50	25405-059070	-	25405-059270	25405-059370	25405-059570
		100	25405-109070	-	25405-109270	25405-109370	25405-109570
		150	25405-159070	-	25405-159270	25405-159370	25405-159570
		250	25405-259070	-	25405-259270	25405-259370	25405-259570


Hypersil GOLD Phenyl

Contains a C4 linker which allows for superior alignment of the phenyl ring with aromatic molecules

- Enhanced pi-pi interactions with aromatics
- Moderate hydrophobicity
- Outstanding peak shape and sensitivity



Antidepressants



Hypersil GOLD Phenyl 1.9µm, 50 x 2.1mm

Mobile Phase A:	0.1% Formic acid
Mobile Phase B:	0.1% Formic acid in MeCN
Gradient:	10 – 60% B in 3.4mins, 60 – 90% B in 0.24 min
Temperature:	60°C
Flow Rate:	0.5mL/min
Injection Volume:	0.7µL
Detection:	UV, 225 and 254nm
Analytes:	1. Uracil 2. Acetaminophen 3. p-Hydroxybenzoic acid 4. o-Hydroxybenzoic acid 5. Oxazepam 6. Diazepam 7. Di-isopropyl phthalate 8. Di-n-propyl phthalate


Hypersil GOLD Phenyl

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	25902-052130	-	-	-
		100	-	25902-102130	-	-	-
		150	-	25902-152130	-	-	-
		200	-	25902-202130	-	-	-
3	Drop-in Guard (4/pk)	10	25903-011001	25903-012101	25903-013001	25903-014001	25903-014001
	HPLC Column	50	-	25903-052130	-	-	-
		100	-	25903-102130	25903-103030	-	25903-104630
		150	25903-151030	-	25903-153030	-	25903-154630
5	Drop-in Guard (4/pk)	10	-	25905-012101	25905-013001	25905-014001	25905-014001
	HPLC Column	50	-	25905-052130	-	-	25905-054630
		100	-	25905-102130	25905-103030	-	-
		150	-	-	-	-	25905-154630
		250	-	-	25905-253030	25905-254030	25905-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00


Hypersil GOLD CN

For both normal phase and reversed-phase separations

- Provide alternative selectivity with lower hydrophobicity
- Excellent peak shape
- Outstanding sensitivity
- Less retention for faster analysis



Organic acids



Hypersil GOLD CN 5µm, 150 x 4.6mm

Mobile Phase A:	25 mM KH ₂ PO ₄ pH2
Mobile Phase B:	MeOH
Temperature:	25°C
Flow Rate:	1.5mL/min
Detection:	UV, 230nm
Analytes:	1. 4-Fluorobenzoic 2. o-Toluic Acid 3. p-Toluic Acid 4. 2,4,6-Trimethylbenzoic Acid 5. 2,5-Dimethylbenzoic Acid


Hypersil GOLD CN

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	25802-052130	-	-	-
		100	-	25802-102130	-	-	-
		200	-	25802-202130	-	-	-
3	Drop-in Guard (4/pk)	10	25803-011001	25803-012101	25803-013001	25803-014001	25803-014001
	HPLC Column	30	-	25803-032130	-	-	-
		50	-	25803-052130	-	-	-
		100	-	25803-102130	25803-103030	-	25803-104630
		150	25803-151030	25803-152130	25803-153030	-	25803-154630
5	Drop-in Guard (4/pk)	10	-	25805-012101	25805-013001	25805-014001	25805-014001
	HPLC Column	50	-	25805-052130	-	-	25805-054630
		100	-	25805-102130	25805-103030	-	-
		150	-	25805-152130	-	-	25805-154630
		250	-	25805-252130	25805-253030	25805-254030	25805-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00


Hypersil GOLD Amino

A high performance aminopropyl phase that gives excellent chromatographic properties in three modes: weak anion exchange, reversed-phase and normal phase

- Retains anions and organic acids in weak anion exchange
- Excellent for carbohydrate analysis in reversed-phase
- Alternative selectivity to silica columns in normal phase chromatography
- Outstanding peak shape and sensitivity



Sugars



Hypersil GOLD Amino 5µm, 150 x 4.6mm

Mobile Phase:	MeCN/Water (80:20)
Temperature:	35°C
Flow Rate:	1.2mL/min
Injection Volume:	20µL
Detection:	RI
Analytes:	1. Fructose 2. Glucose 3. Sucrose 4. Maltose 5. Lactose


Hypersil GOLD Amino

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	25702-052130	-	-	-
		100	-	25702-102130	-	-	-
		150	-	25702-152130	-	-	-
		200	-	25702-202130	-	-	-
3	Drop-in Guard (4/pk)	10	25703-011001	25703-012101	25703-013001	25703-014001	25703-014001
	HPLC Column	30	-	25703-032130	-	-	-
		50	-	25703-052130	-	-	-
		100	-	25703-102130	25703-103030	-	25703-104630
		150	25703-151030	25703-152130	25703-153030	-	25703-154630
5	Drop-in Guard (4/pk)	10	-	25705-012101	25705-013001	25705-014001	25705-014001
	HPLC Column	50	-	25705-052130	-	-	25705-054630
		100	-	25705-102130	25705-103030	-	-
		150	-	-	-	-	25705-154630
		250	-	25705-252130	25705-253030	25705-254030	25705-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00


Hypersil GOLD AX

A novel polymeric amine ligand bonded to highly pure base deactivated silica

- Weak anion exchange phase for multiple charged species
- Suitable for HILIC retention and separation of highly polar molecules
- Higher efficiency than polymer based ion exchange columns
- Outstanding peak shape and selectivity



Vitamin C



Hypersil GOLD AX 5µm, 100 x 4.6mm

Mobile Phase: 100 mM Ammonium Acetate pH 6.8/
MeCN (30:70)
Temperature: 30°C
Flow Rate: 0.5 mL/min
Injection Volume: 50 µL
Detection: UV, 240nm


Hypersil GOLD AX

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	26102-052130	-	-	-
		100	-	26102-102130	-	-	-
		150	-	26102-152130	-	-	-
		200	-	26102-202130	-	-	-
3	Drop-in Guard (4/pk)	10	26103-011001	26103-012101	26103-013001	26103-014001	26103-014001
	HPLC Column	30	-	26103-032130	-	-	-
		50	-	26103-052130	-	-	-
		100	-	26103-102130	-	-	26103-104630
		150	26103-151030	26103-152130	26103-153030	-	26103-154630
5	Drop-in Guard (4/pk)	10	-	26105-012101	26105-013001	26105-014001	26105-014001
	HPLC Column	50	-	-	-	-	26105-054630
		100	-	26105-102130	-	-	26105-104630
		150	-	26105-152130	-	-	26105-154630
		250	-	26105-252130	26105-253030	26105-254030	-
		UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00


Hypersil GOLD SAX

A highly stable quaternary amine strong anion exchange column for aqueous and low pH mobile phases

- High stability to aqueous and low pH mobile phases
- Ideally suited to the analysis of smaller organic molecules including nucleotides and organic acids
- Outstanding peak shape and sensitivity



Monophosphates



Hypersil GOLD SAX 5µm, 150 x 4.6mm

Mobile Phase: Aqueous KH_2PO_4 (50 mM, pH 3)
 Temperature: 40°C
 Flow Rate: 1.0mL/min
 Injection Volume: 10µL
 Detection: UV, 254nm


Hypersil GOLD SAX

Particle Size (µm)	Format	Length (mm)	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	100	26302-102130	-	-	-
		150	26302-152130	-	-	-
3	Drop-in Guard (4/pk)	10	26303-012101	26303-013001	26303-014001	26303-014001
	HPLC Column	50	26303-052130	-	-	-
		100	26303-102130	26303-103030	-	26303-104630
		150	-	26303-153030	-	26303-154630
5	Drop-in Guard (4/pk)	10	26305-012101	26305-013001	26305-014001	26305-014001
	HPLC Column	50	26305-052130	-	-	-
		100	26305-102130	26305-103030	-	26305-104630
		150	26305-152130	-	-	26305-154630
		250	26305-252130	26305-253030	26305-254030	26305-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	852-00	852-00	850-00	850-00


Hypersil GOLD Silica

Unbonded, highly pure base deactivated silica media that is the backbone of the Hypersil GOLD range of columns

- Highly pure base deactivated silica media
- Outstanding peak shape and sensitivity



Steroids



Hypersil GOLD Silica 5µm, 150 x 4.6mm

Mobile Phase:	19:1 (v/v) n-C6H14/EtOH
Temperature:	30°C
Flow Rate:	1.5mL/min
Injection volume:	5µL
Detection:	UV, 254nm
Analytes:	<ol style="list-style-type: none"> 1. Progesterone 2. 21-Hydroxyprogesterone-21-acetate 3. 17-α-Hydroxyprogesterone 4. Cortisone 5. 11-α-Hydroxyprogesterone 6. Corticosterone 7. Hydrocortisone


Hypersil GOLD Silica

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	25102-052130	-	-	-
		100	-	25102-102130	-	-	-
		150	-	25102-152130	-	-	-
		200	-	25102-202130	-	-	-
3	Drop-in Guard (4/pk)	10	25103-011001	25103-012101	25103-013001	25103-014001	25103-014001
	HPLC Column	30	-	25103-032130	-	-	25103-034630
		50	-	25103-052130	-	-	-
		100	-	25103-102130	25103-103030	-	-
		150	25103-151030	25103-152130	25103-153030	-	25103-154630
5	Drop-in Guard (4/pk)	10	-	25105-012101	25105-013001	25105-014001	25105-014001
	HPLC Column	50	-	25105-052130	-	-	-
		100	-	25105-102130	25105-103030	-	25105-104630
		150	-	25105-152130	-	-	25105-154630
		250	-	25105-252130	-	25105-254030	25105-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00


Hypersil GOLD HILIC

Retains and separates polar analytes that are problematic using reversed-phase columns

- Alternative selectivity to C18
- Improved sensitivity for MS detection
- Alternative to ion-pair or derivatization
- Outstanding peak shape and selectivity



Urea



Hypersil GOLD HILIC 5µm, 150 x 4.6mm

Mobile Phase:	H ₂ O/MeCN (10:90) + 0.1% formic acid
Temperature:	30°C
Flow Rate:	0.6mL/min
Injection Volume:	1µL
Detection:	+ESI
Analytes:	1. Urea

Hypersil GOLD HILIC

Particle Size (µm)	Format	Length (mm)	1.0mm ID	2.1mm ID	3.0mm ID	4.0mm ID	4.6mm ID
1.9	UHPLC Column	50	-	26502-052130	-	-	-
		100	-	26502-102130	-	-	-
		150	-	26502-152130	-	-	-
3	Drop-in Guard (4/pk)	10	26503-011001	26503-012101	26503-013001	26503-014001	26503-014001
	HPLC Column	30	-	26503-032130	-	-	-
		50	-	26503-052130	-	-	-
		100	-	26503-102130	26503-103030	-	26503-104630
		150	26503-151030	26503-152130	26503-153030	-	26503-154630
5	Drop-in Guard (4/pk)	10	-	26505-012101	26505-013001	26505-014001	26505-014001
	HPLC Column	50	-	26505-052130	-	-	26505-054630
		100	-	26505-102130	26505-103030	-	-
		150	-	-	-	-	26505-154630
		250	-	26505-252130	26505-253030	26505-254030	26505-254630
	UNIGUARD Drop-in Guard Cartridge Holder	10	851-00	852-00	852-00	850-00	850-00